

OPTIFLEX KNEE CONTINUOUS PASSIVE MOTION MACHINE INSTRUCTIONS


BASIC INSTRUCTIONS

1. Make sure the CPM is plugged into the wall. Make sure that the electrical socket used is not operated by a light switch. If it is, make sure the light switch is in the ON position.
2. Turn the Master Power Switch to the ON position. (Located next to where the power cord enters the base of the CPM.)
3. Ensure that the KCPM is at its flattest position prior to placing the patient's leg into the machine.
4. Apply required Safety Straps around the patient's leg. The Patient Service Representative that delivered the KCPM should have showed you how to apply these straps.
5. Double check the controller settings. If they are at the proper settings, push the red start button located at the bottom of the controller to start the movement.




CHANGING SETTINGS ON THE CONTROLLER

1. The EXTENSION or how flat your leg gets should be set by the Patient Service Representative prior to leaving your home or Facility. Should it need to be adjusted press the button with the picture of a flat leg located on the left side of the controller, while keeping this button depressed, press the up or down buttons to adjust.



CHANGING SETTINGS ON THE CONTROLLER

2. The FLEXION or how much your leg bends can be adjusted by pressing the button that has a picture of a bent leg on the left side of the controller. While keeping this button pressed in, push the up or down button to adjust.



Please remember that use of the CPM does not take the place of your Physical Therapy Exercises. If you any issues with the CPM, or have any questions call 1-800-996-4001.

Doctor Protocols: Start _____ to _____, Increase _____ to _____ per day, End _____