

TITAN™ WRIST - LACING ORTHOSIS

68.25/95.61H/51.19UW PDAC Approved L3908



- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Universal right and left sizes reduce inventory needs
- Adjustable lacing system handles a wide range of wrist sizes
- Palmal and dorsal stays are malleable for customization and removable for "step-down" rehabilitation
- Stockinette tongue eases application
- Anatomically designed for a better fit

Order No. Description / Options

Titan™ Wrist 1365.659.0611 450 Titan™ Wrist & Forearm 1815.1721.3451 452

Specify Right or Left

XS & XL sizes also available for 450 only

TITAN™ WRIST SIZING:

Universal - Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm) **XS** - Wrist Circumference 4.0" - 6.0" (10.2cm - 15.2cm) **XL -** Wrist Circumference 7.5" - 11.0" (19.1cm - 27.9cm)



Dorsal stay pod can be moved proximal or distal for desired level of dorsiflexion restriction





TINY TITAN™ WRIST

 Pediatric splint constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior

 Anatomically contoured palmal stay is malleable for customization and removable for "step-down" rehabilitation

■ Stockinette tongue eases application

Anatomically designed for a better fit

Order No. Description / Options **450 TINY** Tiny Titan™ Wrist

Specify Right or Left 1365.659.0611





Wrist Circumference 3.5" - 5" (8.9cm - 12.7cm)





SUEDE LACING WRIST ORTHOSIS

PDAC Approved L3908 68.25/95.61H/51.19UW

- Constructed of perforated black suede laminated to moisture moving poly-felt
- Contoured design provides secure fit and allows full finger and thumb function
- Lacing design provides for easy secure closures
- Stockinette tongue eases application
- Anatomically designed for a low profile fit
- Aluminum palmal and dorsal stays are malleable and removable

Order No. Description / Options

434 Lacing Wrist Orthosis 1395.779.6811 444 Lacing Wrist & Forearm 1695.7811.1441

Specify Right or Left

SUEDE LACING WRIST SIZING: Wrist Circumference





UNIVERSAL WRIST ORTHOSIS





- Palmal and dorsal malleable aluminum stays may be molded for a custom fit or removed for stepdown rehabilitation
- Universal right and left sizes reduce inventory needs
- Elastic carpal strap provides support & comfort

Order No. Description / Options

438 Universal Wrist Orthosis 1225.858.1401439 Universal Wrist & Forearm 1395.779.6811

Specify Right or Left

UNIVERSAL WRIST SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm)



PHOMFIT™ WRIST ORTHOSIS

PDAC Approved L3908 68.25/95.61H/51.19UW

- Comfortable foam and nylon tricot lining, laminated to a sturdy black polyester outer fabric provides comfort and durability
- Removable aluminum palmal stay is malleable
- Stockinette tongue and dorsal hook & loop closures offer simple application
- No exposed lateral hook

Order No. Description / Options

436 Phomfit™ Wrist Orthosis - Short 1425.899.1121

437 Phomfit[™] Wrist Orthosis - Long

Specify Right or Left

PHOMFIT™ WRIST SIZING: Wrist Circumference

XS: 4.5" - 5.5" (11.4cm - 14.0cm)
Sm: 5.0" - 6.5" (12.7cm - 16.5cm)
Med: 6.0" - 7.5" (15.2cm - 19.1cm)
Lg: 7.0" - 8.5" (17.8cm - 21.6cm)
XL: 8.0" - 9.5" (20.3cm - 24.1cm)





UNO WHO® - WRIST HAND ORTHOSIS

PDAC Approved L3908 68.25/95.61H/51.19UW



- Movable and malleable stay pod provides custom sizing and fit
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability

Order No. Description / Options

428 UNO Wrist Hand Orthosis® 1455.9101.7321

Specify Right or Left

UNO WHO® WRIST SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm)



Movable stay pod









GEL WRIST WRAP

428

- Kuhl[™] perforated neoprene construction provides comfortable compression
- Gel pad cushions and protects the wrist from shock and vibration
- Gel visco-elastic properties move in a "fluid-like" manner with the skin
- Patented patient assist strap makes one-handed application easy
- Ideal for reducing scarring and increasing range of motion following Carpal Tunnel Surgery
- Universal sizing reduces inventory

Order No. Description / Options

3846 Gel Wrist Wrap (Gel included) 2995.7902.6452

GEL Gel Pad Insert 1705.4911.9441

Specify Right or Left

GEL WRIST SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm)





MODABBER™ WRIST ORTHOSIS

PDAC Approved L3908 68.25/95.61H/51.19UW



- Patented patient assist single strap for easy application
- Kuhl[™] perforated neoprene construction provides comfortable compression
- 4 way stretch neoprene ensures brace conforms to contour of hand
- Aluminum palmal stay is removable and malleable for a custom fit
- Universal sizing reduces inventory needs

Description / Options Order No.

Kuhl™ Modabber™ Wrist - Long 1695.7811.1441 5818 $\mathsf{Kuhl}^\mathsf{TM}\ \mathsf{Modabber}^\mathsf{TM}\ \mathsf{Wrist}\ \text{-}\ \mathsf{Short}\ \ 1495.7401.1721$ 5819

5818 PED Kuhl[™] Modabber[™] Wrist - Long Pediatric **5819 PED** Kuhl[™] Modabber[™] Wrist - Short Pediatric

Specify Right or Left

MODABBER™ WRIST SIZING:

Wrist circumference

Pediatric 4.0" - 6.0" (10.2cm - 15.2cm) Standard 4.5" - 9.5" (11.4cm - 24.1cm)





M. Ramin Modabber, M.D.

Dr. Modabber is a Hand & Upper Extremity Specialist at Santa Monica Orthopaedic and Sports Medicine Group and Saint John's Health Center in Santa Monica, CA



DYNADIGIT®

PDAC Approved L3807 207.39/214.77UW



- Neoprene generates dynamic action at any joint
- Apply hook dorsally to restrict MCP motion
- Interchangeable stays fit in dorsal pocket
- Malleable stay immobilizes
- Coiled stay assists extension or resists flexion
- Easily attaches to the Modabber neoprene wrist cuff for adjustable tension

Available individually or with the Modabber brace for wrist control.

Order No. Description 2695.7881.1922 Modabber™ Wrist with DynaDigit® 5821

Specify Right or LeftDynaDigit® - Single Unit 605A

1280.698.8801 MODABBER™ WRIST SIZING:

Universal: Wrist circumference 4.5" - 9.5" (11.4 cm - 24.1cm) DYNADIGIT® SIZING:

Universal - Fits most finger lengths



US Patents 6,142,966

Michael Behrman, M.D.

Developed in collaboration with Michael Behrman, M.D. of Associated Hand Surgeons, Santa Barbara, CA.



TITAN™ THUMB ORTHOSIS

PDAC Approved L3807 207.39/214.77UW



- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Adjustable lacing pull easily fits a wide range of wrist sizes
- Malleable oblique thumb stay allows desired thumb position and stability
- Provides abduction control, angle adjustability, and easy reapplication
- Movable stay pod makes this brace universal and reduces inventory
- Stockinette tongue and dorsal loop & lock closures offer simple application
- The oblique stay, the short malleable dorsal stay and long radial stay offer cone like compression

Order No. Description / Options

455 Titan[™] Thumb Orthosis 2135.5941.5181 Long Titan[™] Thumb Orthosis 2565.6971.0812

Specify Right or Left

XS & XL sizes also available for 455 only

TITAN™ THUMB SIZING:

Universal - Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm)

XS - Wrist Circumference 4.0" - 6.0" (10.2cm - 15.2cm)

XL - Wrist Circumference 7.5" - 11.0" (19.1cm - 27.9cm)



The oblique thumb stay can be bent to allow desired thumb positioning

TINY TITAN™ THUMB





Anatomically designed for a better fit

Order No. Description
455 TINY Tiny Titan™ Wrist

Specify Right or Left 2135.5941.5181



): 9cm - 12 7cm)





Wrist Circumference 3.5" - 5" (8.9cm - 12.7cm)





UNO WHT® - WRIST HAND THUMB

PDAC Approved L3807 207.39/214.77UW



- Movable stay pod makes this brace universal reducing inventory
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Useful in treatment of de Quervain's tendinitis, Gamekeeper injuries, Scaphoid/Navicular injuries, and is used for early cast removal
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability
- Malleable oblique thumb stay allows desired thumb position and stability
- The oblique stay, the short malleable dorsal stay and long radial stay offer cone like compression

Order No. Description / Options

429 UNO Wrist Hand Thumb™ Orthosis


Specify Right or Left 2295.7061.1591

UNO WHT® WRIST SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm)

Removable, malleable Aluminum Palmal Stay











PHOMFIT™ WRIST HAND THUMB ORTHOSIS

PDAC Approved L3807 207.39/214.77UW

- Aluminum removable palmal stay allows for custom fit and adjustment
- Stockinette tongue and dorsal hook & loop closures offer simple application
- Radial thumb stay is malleable
- Useful in treatment of de Quervain's Tendinitis, Gamekeeper Injuries, Scaphoid/Navicular Injuries, and is used for early cast removal
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability

Order No. Description / Options

447 Phomfit™ Wrist Hand Thumb Orthosis Specify Right or Left 2395.7761.6302



 XS:
 4.5" - 5.5" (11.4cm - 14.0cm)

 Sm:
 5.0" - 6.5" (12.7cm - 16.5cm)

 Med:
 6.0" - 7.5" (15.2cm - 19.1cm)

 Lg:
 7.0" - 8.5" (17.8cm - 21.6cm)

 XL:
 8.0" - 9.5" (20.3cm - 24.1cm)







SANTA BARBARA THUMB SPLINT

PDAC Approved L3923 74.31



- Radial, palmal and dorsal pockets contain malleable aluminum and spiral stays provide customized support
- Ideal for treating Gamekeeper's/Skier's Thumb, Basal Joint Arthritis, and Thumb Tendonitis
- Wrap around design for one-handed application
- Fits both the right and left hand

Order No. Description

3842 Santa Barbara Thumb Splint 1910.7331.4261

SANTA BARBARA THUMB SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm) Thumb Hole Sizing 2.5" - 3.375" (6.4cm - 8.6cm)





Michael Behrman, M.D.

Developed in collaboration with Michael Behrman, M.D. of Associated Hand Surgeons, Santa Barbara, CA. He is a fellow of the American Academy of Orthopedic Surgeons and member of the American Society for Surgery of the Hand.



TRIMABLE™ THUMB ORTHOSIS

- 1/8" (3mm) Nylon covered neoprene provides heat and compression
- Melco taping in thumb section for trimming
- Wrap around design with built in assist strap for easy application
- Fits both the right and left hands

Order No. Description

Trimable™ Thumb Orthosis 1895.7231.1161 3810



TRIMABLE™ THUMB SIZING:

Thumb Hole Sizing 2.125" - 2.375" (5.4cm - 6.0cm) REG: 2.375" - 3.0" (6.0cm - 7.6cm)

3.0" - 3.625" (7.6cm - 9.2cm)

Wrist Circ. 4.5"- 5.5" (11.4cm - 14.0cm) 5.0"- 9.5" (12.7cm - 24.1cm) 5.0"- 9.5" (12.7cm - 24.1cm)





NEOPRENE THUMB ORTHOSIS

PDAC Approved L3923 74.31

- Universal size one size fits all
- Neoprene construction provides an excellent fit of the hand and thumb
- Elasticized straps for comfortable compression
- Supports the MP and CMC joints while allowing full finger motion
- An excellent resting splint for thumb tendonitis and de Quervain's tendonitis

Order No. Description

Neoprene Thumb Orthosis 1425.899.1121 3813

NEOPRENE THUMB SIZING:

Wrist Circumference 4.5" - 9.5" (11.4cm - 24.1cm) Thumb Hole Sizing 2.5" - 3.375" (6.4cm - 8.6cm)





MODABBER™ THUMB ORTHOSIS

PDAC Approved L3923 74.31



■ Constructed of Kuhl[™] perforated neoprene for cool, comfortable compression

■ Brace fits left or right hand to reduce inventory

Aluminum stay and steel spiral stay are interchangeable and removable

Ideal for treating Basal Joint Arthritis, Throttle Thumb, and Thumb Tendonitis

Order No. Description / Options

5803 $\mathsf{Kuhl}^\mathsf{TM}$ $\mathsf{Modabber}^\mathsf{TM}$ Thumb $\mathsf{Orthosis}$

1910.7331.4261

MODABBER™ THUMB SIZING:

XS:

XL:

Thumb Hole Sizing

2.125" - 2.375" (5.4cm - 6.0cm) 2.375" - 3.0" (6.0cm - 7.6cm) 3.0" - 3.625" (7.6cm - 9.2cm) REG:

Wrist Circ.

4.5"- 5.5" (11.4cm - 14.0cm) 5.0"- 9.5" (12.7cm - 24.1cm) 5.0" - 9.5" (12.7cm - 24.1cm)



Pocket flap makes it possible to interchange and remove stays



Taping above stay allows trimming for freedom of the IP joint



M. Ramin Modabber, M.D.

Dr. Modabber is a Hand & Upper Extremity Specialist at Santa Monica Orthopaedic and Sports Medicine Group and Saint John's Health Center in Santa Monica, CA



CMC CONTROLLER™

- Places the thumb in the optimum position by simultaneously abducting the thumb and distracting (seating) the CMC joint
- Low profile allows normal motion of the thumb while supporting the CMC joint
- Two sizes minimizes inventory (no rights or lefts)

Order No. Description / Options

CMC Controller™ 1495.7401.1721 3800

3801 CMC Controller™ with pressure dot 1895.7231.1161

CMC CONTROLLER SIZING:

XS/S Fits wrist circumference 4.5" - 6.5" (11.4cm - 16.5cm) M/L Fits wrist circumference 6.5" - 10" (16.5cm - 25.4cm)

X-Large size also available



David M. Auerbach, M.D.

Dr. David Auerbach practices at Southern California Orthopedic Institute (SCOI) in Van Nuys, CA.



Shane Woods, COTA/L

Mr. Shane M. Woods is the Directer of Occupation and Hand Therapy for the Hand Center of Ventura County, Oxnard, CA.



TKO® - THE KNUCKLE ORTHOSIS

PDAC Approved L3807 207.39/214.77UW



- Provides positioning for fractures or injuries to metacarpals, phalanges,
 MCP and IP joints
- Universal size, Right and Left
- Straps are infinitely adjustable for fit and comfort (pain & swelling)
- Distal buddy strapping eliminates rotation of the finger
- Padded stays are malleable
- Reversible one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle & Ring Finger splint

Order No. Description / Options **3848** TKO® (The Knuckle Orthosis)

Specify Right or Left 2565.6971.0812

TKO® SIZING:

Measure circumference at palmar crease Regular fits less than 10.5in circumference

X-Large size also available, fits greater than 10.5in circumference XL2775.3491.9532



David Garelick, M.D.

Developed in collaboration with Dr. David Garelick, Sports Medicine Physician at Midwest Orthopaedics Illinois Bone & Joint, Chicago. IL



Edward Cotton, OTC

Mr. Cotton practices at Cook County Hospital Chicago, IL



"SQUEEZE" ULNAR COMPRESSION WRAP

- Compressive and comfortable wrist support effective for Treatment of Ulnar Wrist Pain where compression of the wrist (squeeze) provides relief.
- Universal; fits both right and left wrist
- Easy application
- Non-slip interior prevents slippage without a tourniquet effect
- Perforated to let moisture escape

Order No. Description

2801 Squeeze Ulnar Compression Wrap 1995.7931.6961



Jamie Harrison, OTR, CHT

Developed in collaboration with Ms. Jamie Harrison who is a Certified Hand Therapist at Ortho Montana Orthopedics and Sports Medicine in Billings, MT.




For best result Pressure Dot is applied over the Ulna in supination





UNIVERSAL MATT STRAP™



- Foamed silicone pad absorbs shock
- Trimable strap allows size customization between 7" to 20" (17.8cm to 50.8cm) circumference
- One size fits all

Order No. Description

3709 Universal Matt Strap[™] 1595.7111.6531 3709






UNIVERSAL MATT STRAP™ SIZING

(Encircle the area of application to measure.)

Fits circumference of 7" - 20" (17.8cm - 50.8cm). Simply trim excess material.



$\textbf{MATT STRAP}^{\scriptscriptstyle\mathsf{TM}}$

US Patent D500,137

- Useful for tennis and golfer's elbow
- Foamed silicone pad absorbs shock

Order No. Description

Matt Strap™ 1295.709.1011 3704



Loop & lock design provides focused compression over the shock absorbing buttress



Matthew Shapiro, M.D.

Developed in collaboration with Dr. Matthew Shapiro who specializes in knee, shoulder and elbow injuries out of Slocum Orthopedics in Eugene, OR.



Measure widest part of forearm

XS: 7"-8" (17.78cm - 20.32cm) **Med**: 9"-10.5" (22.86cm - 26.67cm) **Lg**: 10.5"-12" (26.67cm - 30.48cm)

Sm: 8"-9" (20.32cm - 22.86cm) 12"- 14" (30.48cm - 35.56cm) **XXL**: 13"- 16" (33.02cm - 40.64cm)



SEMI-UNIVERSAL MATT STRAP™

- Useful for tennis and golfer's elbow
- Non-loop lock construction
- Foamed silicone pad absorbs shock
- Two sizes

Order No. Description

 $Semi\text{-}Universal\ Matt\ Strap^{^{\text{\tiny TM}}}$ 3707 1495.7401.1721





The elastic double-pull strap provides unique additional compression and adjustment

SEMI-UNIVERSAL MATT STRAP $^{\text{\tiny TM}}$

Measure widest part of forearm

S/M: 11"-14" (27.9cm - 35.6cm) **L/XL**: 14"-16" (35.6cm - 40.6cm)



TENNIS ELBOW - WITH PRESSURE PAD

- Low profile design
- 1/8" (3mm) Neoprene with 2-sided nylon outer covering
- Pressure pad & elastic strap provide compression
- Provides relief for epicondylitis

Order No. Description

3838 Tennis Elbow - Pressure Pad

1195.738.6101

TENNIS ELBOW SIZING:

Measure the widest part of the forearm.

 XS:
 8"-9" (20.3cm - 22.9cm)
 Sm:
 9"-10" (22.9cm - 25.4cm)

 Med:
 10" -11" (25.4cm - 27.9cm)
 Lg:
 11"-12" (27.9cm - 30.5cm)

 XL:
 12"-14" (30.5cm - 35.6cm)
 XXL:
 14"-16" (35.6cm - 40.6cm)



FREEHAND DEX BRACE®

US Patent Pending

PDAC Approved L3762 89.18/92.36UW

- A comfortable, dynamically extending alternative for treatment of Cubital Tunnel Syndrome.
- Dynamic extension returns the flexed elbow to the extended position even during sleep
- Brace is adjustable, flexible, lightweight, low profile with minimal limb contact

Order No. Description

469 Freehand DEX Brace® 3995.7972.6933

FREEHAND DEX SIZING:

S/M Fits most females and smaller male arms

Lg/XL Fits most males and larger female arms



Roy A. Meals, M.D.

Developed in collaboration with Dr. Meals, Orthopedic Hand Surgeon, Inventor of Freehand Braces and a Clinical Professor at UCLA



CUBITAL COMFORT™ BRACE

PDAC Approved L3762 89.18/92.36UW



- Ideal for treatment of Cubital Tunnel Syndrome
- Two aluminum malleable stays positioned inside the elbow prevent excessive elbow flexion
- Aluminum stays are pre-bent at 20 degrees of full flexion
- lacktriangle Two stockinette arm bands aid in arm positioning

Order No. Description

3847 Cubital Comfort™ Brace 2595.7181.6022

US Patents: 7,074,202; D523,961

CUBITAL COMFORT™ BRACE SIZING:

Measure length from mid forearm to mid humerus

Short (10 3/4") Long (13 1/4") Length Less than 11 1/2" 11 1/2" and above



David M. Auerbach, M.D.

Dr. David Auerbach practices at Southern California Orthopedic Institute (SCOI) in Van Nuys, CA.





SHOULDER IMMOBILIZER

128.81/163.73H/96.60UW PDAC Approved L3670

- Padded pressure sensitive straps are adjustable and comfortable
- Body strap is easily adjusted or can be removed
- Deep pocket design
- Two adjustment points for each strap allows for better fit

Order No. Description

501 Shoulder Immobilizer Sling

1280.698.8801

SHOULDER IMMOBILIZER SIZING:

Measure from elbow to metacarpals **XS**: 7.5" (Ht. 19.1cm) x 11"(Lg. 27.9cm) **M**: 7.5" (Ht. 19.1cm) x 15.5" (Lg. 39.4cm) **L**:

XL: 7.5" (Ht. 19.1cm) x 20" (Lg. 50.8cm)

SM: 7.5" (Ht. 19.1cm) x 14" (Lg. 35.6cm) 7.5" (Ht. 19.1cm) x 18" (Lg. 45.7cm)



ARM SLING

- Padded strap offers customized fitting and adjustment
- Two adjustment points allow for more comfortable fit
- Deep pocket design accommodates the largest of casts and dressings
- Open end with thumb loop adds comfort and eliminates slippage

Order No. Description

500 Arm Sling 865.606.537

ARM SLING SIZING:

Measure from elbow to metacarpals **XS**: 7.5" Ht. (19.1cm) x 11" Lg. (27.9cm) **SM**: 7.5" Ht. (19.1cm) x 13.5" Lg. (34.3cm) 7.5" Ht. (19.1cm) x 15.5" Lg. (39.4cm) Ŀ 7.5" Ht. (19.1cm) x 17.5" Lg. (44.5cm) **XL**: 7.5" Ht. (19.1cm) x 20" Lg. (50.8cm)

Order No. Description / Options Universal Arm Sling 506

UNIVERSAL ARM SLING SIZING:

One size fits all



US Patent: 7,189,213;

7,300,410



KNAPP SAK II® - SHOULDER ABDUCTION KIT PDAC Approved L3670 128.81/163.73H/96.60UW

- Indicated for use after Bankhart repairs, glenohumeral dislocation, and subluxation, soft tissue strains and repairs, and immobilization after rotator cuff repairs
- Maintains shoulder in neutral position allowing approximately 30 degrees of shoulder abduction
- Prevents excessive internal rotation of the shoulder
- Relieves tension from the repaired rotator cuff tendon (primarily the supraspinatus)
- Color coded for easy application

Order No. Description 510 Knapp Sak II® 5345.2473.3454

KNAPP SAK II® SIZING: One size fits all



Thomas P. Knapp, M.D.

Developed in collaboration with Dr. Knapp, Santa Monica Orthopaedic and Sports Medicine Group, Santa Monica, CA



SHOULDER HOULDER II

PDAC Approved L3670 128.81/163.73H/96.60UW

- Our breathable sling is constructed of our NoSwet® Mesh material, combined with a soft moisture guick-wicking liner to provide your patient with comfort and support
- Four sizes of 15° Abduction Pillows allow for ambidextrous and custom patient fitting
- The unique body and shoulder straps provide quick adjustment, easy application, and padded comfort
- Quick release buckles allow for easy application and removal of the brace
- A posterior gusset comfortably holds the elbow in the desired position without providing unnecessary pressure to the area

Order No. Description

514

Shoulder Houlder II 5795.7504.6294

SHOULDER HOULDER II SIZING:

Measure from Olecranon to Knuckles:

11.5" - 13" (29cm - 33cm) **M** 13.5" - 14.5" (33cm - 37cm) 15" - 16" (37cm - 41cm) **XL** 16.5" - 17.5" (42cm - 45cm)



SHOULDER HOULDER II NEUTRAL

PDAC Approved L3670 128.81/163.73H/96.60UW

- Our breathable sling is constructed of our NoSwet® Mesh material, combined with a soft moisture quick-wicking liner to provide your patient with comfort and support
- Two sizes of Neutral Pillows allow for ambidextrous and custom patient fitting
- The unique body and shoulder straps provide guick adjustment, easy application, and padded comfort
- Quick release buckles allow for easy application and removal of the brace
- A posterior gusset comfortably holds the elbow in the desired position without providing unnecessary pressure to the area

Order No. Description

Shoulder Houlder II Neutral 524

5795.7504.6294

SHOULDER HOULDER II NEUTRAL SIZING:

Measure from Olecranon to Knuckles:

11.5" - 13" (29cm - 33cm) **M** 13.5" - 14.5" (33cm - 37cm)

15" - 16" (37cm - 41cm) **XL** 16.5" - 17.5" (42cm - 45cm)



VAPOR IMMOBILIZATION SLING

PDAC Approved L3670 128.81/163.73H/96.60UW

- Our breathable sling is constructed of our NoSwet® Mesh material, combined with a soft moisture quick-wicking liner to provide your patient with comfort and support
- The unique body and shoulder straps provide quick adjustment, easy application, and padded comfort
- Quick release buckles allow for easy application and removal of the brace
- A posterior gusset comfortably holds the elbow in the desired position without providing unnecessary pressure to the area
- Color coded strap/buckle for easier application
- Posterior O-ring on sling accommodates patient allowing for better fit and comfort

Description

504 Vapor Immobilization Sling

2995.7902.6452

VAPOR IMMOBILIZATION SLING SIZING:

Measure from Olecranon to Knuckles:

11.5" - 13" (29cm - 33cm) **M** 13.5" - 14.5" (33cm - 37cm) 15" - 16" (37cm - 41cm) **XL** 16.5" - 17.5" (42cm - 45cm)





CLAVICLE SPLINT

PDAC Approved L3650 67.55/80.29H/50.67UW

- Ideal for clavicle fracture splinting and posture control
- Comfortable foam construction provides proper immobilization
- Contoured urethane foam is hook & loop sensitive allowing maximum adjustability
- Two buckles on rear felt pad are used for tension adjustment

Order No. Description

661

Clavicle Splint 1145.208.379

CLAVICLE SPLINT SIZING:

Measure from center of shoulder blades around shoulder:

	Min	Max	Min	Max
XS	15"	21"	38.10cm	53.34cm
SM	19"	28"	48.26cm	71.12cm
M	24"	34"	60.96cm	86.36cm
L	30"	38"	76.20cm	96.52cm
XL	34'	40"	86 36cm	101 60cm



CERVICAL COLLAR

30.89/45.83H/23.17UW PDAC Approved L0120



- Medium/firm density urethane foam provides excellent support
- Comfortable, breathable natural stockinette cover
- Vinyl support strap with extended contact closure for a universal fit and classic look
- Available in 2", 3", and 4" heights

Order No. Description

2" Cervical Collar 955.966.218 619

3" Cervical Collar 629 4" Cervical Collar 639



TWO-PULL 'LASTIC LUMBAR SUPPORT

- Manufactured of heavy elastic material
- Two side pulls provide added support
- 4 bone stays provide protection against brace rollover
- Pulls contain assist loops to aid in patient application

Order No. Description

151 Two-Pull 'Lastic Lumbar Support

2675.3781.4722 2895.7202.1642XXL



TWO-PULL 'LASTIC LUMBAR SUPPORT SIZING:

Hip Sizing for 151:

25" - 29" (63.5cm - 73.7cm) XS SM 29" - 33" (73.7cm - 83.8cm) 33" - 37" (83.8cm - 94.0cm) 37" - 41" (94.0cm - 104.1cm) L 41" - 45" (104.1cm - 114.3cm) XXL 45" - 49" (114.3cm - 124.5cm)



LEI BACK™ II

374.59/388.76UW PDAC Approved L0627

- 1 unique Back Orthosis design offered in 2 distinct materials for patient support & comfort choices
- Kuhl™ perforated neoprene provides comfortable compression and support
- No-Swet® spacer fabric provide breathability, comfort and support
- Posterior polypro plastic panel contours to personalized anatomical fit and provides rigid
- A posterior secret pocket is provided for the lordotic foam pad which offers additional comfort, and step-down support opportunities
- Anterior polypro plastic panel provides additional inner abdominal pressure and support while offering step-down capability when less support is appropriate
- Presented in an exceptional personal brace carrying case providing value and utility Order No. Description

Kuhl[™] Lei Back[™] II 9995.7996.<u>69</u>48 5805

7805 No Swet® Lei Back™ II

LEI BACK™ II SIZING:

Measure widest part of hip/waist: 24" - 34" (61.0cm - 86.4cm) 34" - 46" (86.4cm - 116.8cm) 45" - 58" (114.3cm - 147.3cm) XS: **REG**:

XL:





LEI BACK™ LUMBOSACRAL BRACE

PDAC Approved L0625 50.19/52.11UW

- Kuhl[™] perforated neoprene construction provides comfortable compression
- 4-way stretch of neoprene easily conforms to a variety of body contours
- Double pull provides strong support and compression ability
- Universal sizing reduces stocking needs
- 4 flexible spiral stays prevent brace rollover
- Does not restrict motion and generates inner abdominal pressure

Order No. Description

5809

10" (25.4cm) Kuhl™ Lei Back™ Brace 5809

4350.5403.8963

4950.5643.8024XL

5550 5883 8174XXL LEI BACK SIZING:

Measure widest part of hip/waist: 24" - 34" (61.0cm - 86.4cm) 34" - 46" (86.4cm - 116.8cm) XS:

REG: 45"- 58" (114.3cm - 147.3cm) XL:







Plastic Insert



DOUBLE BACK LSO

PDAC Approved L0626 71.02/73.70UW

- Heavy elastic construction ensures comfort
- Elastic pocket accommodates thermoplastic insert with assist strap for easy application
- Double side pulls across front for added support and compression
- Supports the lower back and encourages normal range of motion
- Unique polypro insert provides support

Order No. Description / Options

Double Back LSO with Posterior Plastic Insert 3995.7972.6933 103

DOUBLE BACK SUPPORT SIZING:

Measure widest part of hip/waist:

XS: 32" - 35" (81.3cm - 88.9cm) **Sm**: 36" - 41" (91.4cm - 104.1cm) **Med**: 42" - 45" (106.7 cm - 114.3 cm) **Lg**: 46" - 50" (116.8 cm - 127.0 cm) **XX**: 51" - 55" (129.5 cm - 139.7 cm) **XXL**: 56" - 60" (142.2 cm - 152.4 cm)







3709

3704

UNIVERSAL MATT STRAP™



- Trimable strap allows size customization between 7" to 20" (17.8cm to 50.8cm) circumference
- Shock absorption buttress
- One size fits all



Daniel .

Sm:

Lg:


XXL:



UNIVERSAL MATT STRAP™ SIZING

(Encircle the area of application to measure.) Fits circumference of 7" - 20" (17.8cm - 50.8cm). Simply trim excess material.





US Patent D500,137

Works for arm, knee, upper arm and upper leg tendonitis

- Useful for patella tendonitis, "Jumpers Knee", and minor knee pain
- Shock absorption buttress

Order No. Description

Matt Strap[™] 1295.709.1011 3704



for sizing and custom fitting

Matthew Shapiro, M.D.

Dr. Shapiro specializes in knee, shoulder and elbow iniuries.



Loop & lock design provides focused compression over the shock absorbing buttress

MATT STRAP™

Measure circumference of affected area

7" - 8" (17.78cm - 20.32cm) **Med**: 9" - 10.5" (22.86cm - 26.67cm)

12" - 14" (30.48cm - 35.56cm)

8" - 9" (20.32cm - 22.86cm) 10.5" - 12" (26.67cm - 30.48cm)

13" - 16" (33.02cm - 40.64cm)





The elastic double-pull strap provides unique additional compression and adjustment

3707

SEMI-UNIVERSAL MATT STRAP™

- Works for arm, knee, upper arm and upper leg tendonitis
- Non-loop lock construction
- Unique double pull compression straps
- Two sizes

Order No. Description

Semi-Universal Matt Strap™ 1495.7401.1721 3707





The elastic double-pull strap provides unique additional compression and adjustment

SEMI-UNIVERSAL MATT STRAP™ SIZING

Measure circumference of affected area

S/M: 11" - 14" (27.9cm - 35.6cm) **L/XL**: 14" - 16" (35.6cm - 40.6cm)



SHIELDS™ BRACE

- Kuhl[™] perforated neoprene construction provides comfortable compression
- Dynamic patella stabilization provided through entire range of motion



US Patent: 5,759,167

- Independent calf and thigh strap application adjusts for most leg contours
- Low profile 8 1/2" (21.6cm) design prevents brace slippage
- Successfully controls malalignments of glide, medial/lateral tilt, and anterior/posterior tilt
- Easy-to-apply anterior closure
- The Shields Brace will control the patella during athletic competition

Order No. Description / Options Kuhl™ Shields™ Brace 5345.2473.3454 5674

SHIELDS™ KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 10" - 12" (25.4cm - 30.5cm)
 Sm: 12" - 14" (30.5cm - 35.6cm)

 Med:
 14" - 16" (35.6cm - 40.6cm)
 Lg: 16" - 18" (40.6cm - 45.7cm)

 XL:
 18" - 20" (45.7cm - 50.8cm)
 XXL: 20" - 24" (50.8cm - 61.0cm)

 Larger and custom sizes available upon request.



Developed in collaboration with Clarence L. Shields, MD, Kerlan-Jobe Orthopaedic Clinic, Los Angeles, CA





Neoprene Condyle Pads

comfort and reduce chaffing.

provide protection and

The New Shields Buttress The Trakaderm™ Buttress is

constructed of Polyurethane and is an adherent, soft, gel like material. It helps assist in the active relocation of the patella and places the bone in its pain free position.

HINGED SHIELDS™ BRACE

PDAC Approved L1820 138.11/147.56H/123.80UW

- 1/8" Kuhl[™] perforated neoprene construction for cool compression
- Dynamic patella stabilization provided through entire range of motion
- Independent calf and thigh strap application adjusts for most leg contours
- Successfully controls malalignments of glide, medial/lateral tilt, and anterior/posterior tilt
- Controls patella during strenuous exercise or athletic activity
- The buttress is independent of the support wrap providing dynamic patella stabilization
- Directional Force Flap holds the Buttress in place against the patella
- Neoprene Condyle Pads provide protection and comfort and reduce chaffing
- Wide bi-axial hinges enhance medial/lateral stability and prevent hyperextension provide protection, comfort and reduce chaffing.

Order No. Description / Options **5675** Kuhl™ Hinged Shields™ Brace

7880.6155.8966

SHIELDS™ KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 10" - 12" (25.4cm - 30.5cm)
 Sm: 12" - 14" (30.5cm - 35.6cm)

 Med:
 14" - 16" (35.6cm - 40.6cm)
 Lg:
 16" - 18" (40.6cm - 45.7cm)

 XL:
 18" - 20" (45.7cm - 50.8cm)
 XXL:
 20" - 24" (50.8cm - 61.0cm)

Larger and custom sizes available upon request.









Range of Motion hinges



Range-of-Motion Options:

Extension Stops allow 10°, 20°, 30°, 40° and 50° settings, (preset @ 10° upon delivery)

Flexion Stops at 50°, 60°, 70°, 80°, 90° and 100° increments

Brace can be locked at full extension

application video @hely-weber.com

VELOCITY PS

633.79/662.17H/530.77UW **PDAC Approved L1832**

- Superior knee brace that provides the control and support your patients need while maintaining the flexibility and the restriction of movement that your rehabilitation protocol requires
- Constructed in Kuhl[™] perforated neoprene which breathes well and offers warmth and maximum compression
- Range of Motion Hinges feature a polycentric adjustable design to control extension and flexion with comfortable shock absorbing condylar pads and shells which ensure patient comfort and compliance
- Circumferential straps have elastic hamstring tension relief allowing no interruption in the firing mechanism of the quadriceps muscles while still providing low profile support.
- The brace transforms from left to right by simply moving the patella force strap from the left to right condylar location
- The brace offers effective support for ligament / meniscus tears while having the patented "U-shaped" patella locator buttress and adjustable patella control of the Shields[™] design

Order No. Description / Options

5645PS 12" Velocity Knee, Range of Motion hinges with Patella Stabilization - Pull Up

11995.7938.69101



Comfortable patented "U-shaped" natella locator buttress

Hook side of buttress





VELOCITY PS SIZING:

Measure circumference 6" above mid-patella with leg fully extended.

13" - 16" (33.0cm - 40.6cm) **Sm**: 15" - 18" (38.1cm - 45.7cm) **Lg**: 18" - 20" (45.7cm - 50.8cm) **XXL**: 21" - 24" (53.3cm - 61.0cm) Med: 16" - 19" (40.6cm - 48.3cm) 19" - 23" (48.3cm - 58.4cm) **XXXL**: 22" - 25" (55.9cm - 63.5cm)



PATELLA STABILIZER - Universal Tubular Buttress

- 12" (30.5cm) Nylon covered neoprene base offers comfort, heat and compression
- Medial and lateral spiral stays limit bunching, while allowing full range of motion
- Universal tubular buttress is movable/removable for individualized patient care

Description / Options 3845.2962.8623

3630 1/8" (3mm) Stabilizer - Universal Tubular Buttress, Open Popliteal 3730 1/8" (3mm) Stabilizer - Universal Tubular Buttress, Closed Popliteal

SLIP-ON KNEE SIZING:

Measure at mid-patella with leg fully extended.

12" - 13" (30.5cm - 33.0cm) **Sm**: 13"-14" (33.0cm - 35.6cm) **Med**: 14" - 15" (35.6cm - 38.1cm) **Lg**: 15"-17" (38.1cm - 43.1cm) **XXL**: 18"-20" (45.7cm - 50.8cm) 16" - 18" (40.6cm - 45.7cm)

Larger and custom sizes available upon request.

- 3630 base with an easy to apply **Posterior Closure** design
- Constructed of Kuhl[™] perforated neoprene for cool, comfortable compression

Order No. Description / Options 4095.7682.1843 **3630PC** Posterior Closure Stabilizer - Universal Tubular Buttress

POSTERIOR CLOSURE KNEE SIZING:

Measure at mid-patella with leg fully extended.

XS: 10" - 13" (25.4cm - 33.0cm) **Med**: 14" - 17" (35.6cm - 43.1cm) **Sm**: 12"-15" (30.5cm - 38.1cm) 16"-19" (40.6cm - 48.3cm) Lg: 18" - 22" (45.7cm - 55.9cm) **XXL**: 20"-24" (50.8cm - 61.0cm)

Larger and custom sizes available upon request.









Inferior tubular buttress provides patella control

3641 / 3741

PATELLA STABILIZER - INFERIOR U TUBULAR BUTTRESS

- 12" (30.5cm) Nylon covered neoprene base offers comfort, heat and compression
- Medial and lateral spiral stays limit bunching, while allowing full range of motion
- Sewn in Inferior U Tubular Butress

Order No. Description / Options 4035.5282.0343 **3641** Stabilizer - Inferior U Tubular Buttress, Open Popliteal **3741** Stabilizer - Inferior U Tubular Buttress, Closed Popliteal

KNEE SIZING:

Measure at mid-patella with leg fully extended.

Larger and custom sizes are available upon request.



AXIS HINGED KNEE SLEEVE

PDAC Approved L1810 113.43/126.68H/88.14UW

- 1/8" (3mm) Neoprene with nylon outer covering
- Removable hinges increase treatment options
- Hinges feature hyperextension stops
- Covered hinges allow full sports participation
- Pull up loops aid in application

Order No. Description / Options 4295.7003.1563
3655 12" (30.5cm) Axis Hinged Knee Sleeve - Open Popliteal
12" (30.5cm) Axis Hinged Knee Sleeve - Closed Popliteal

KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 12" - 13" (30.5cm - 33.0cm)
 Sm:
 13" - 14" (33.0cm - 35.6cm)

 Med:
 14" - 15" (35.6cm - 38.1cm)
 Lg:
 15" - 17" (38.1cm - 43.1cm)

 XL:
 16" - 18" (40.6cm - 45.7cm)
 XXL:
 18" - 20" (45.7cm - 50.8cm)

Larger and custom sizes are available upon request.



5694H / 5794H 5690H / 5790H

HINGED LATERAL J PATELLA STABILIZER PDAC Approved L1820 138.11/147.56H/123.80UW

- J Buttress in conjunction with Anterior Flap gives dynamic patella control
- Biaxial Hinges provide medial and lateral stability
- Kuhl[™] perforated neoprene allows moisture to escape and the skin to remain dry

Order No. Description / Options 4650.5523.3593

5694H Kuhl[™] Hinged Lateral J Stabilizer - Open Popliteal with Condyle Pads Kuhl[™] Hinged Lateral J Stabilizer - Closed Popliteal with Condyle Pads

PDAC Approved L1810 113.43/126.68H/88.14UW

Order No.

Description / Options 4935.5543.5914

Kuhl™ Hinged Lateral J Stabilizer - Open Popliteal

Kuhl™ Hinged Lateral J Stabilizer - Closed Popliteal

HINGED LATERAL J PATELLA STABILIZER SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 12" - 13" (30.5cm - 33.0cm)
 Sm:
 13" - 14" (33.0cm - 35.6cm)

 Med:
 14" - 15" (35.6cm - 38.1cm)
 Lg:
 15" - 17" (38.1cm - 43.1cm)

 XL:
 16" - 18" (40.6cm - 45.7cm)
 XXL:
 18" - 20" (45.7cm - 50.8cm)

Larger and custom sizes are available upon request.



KNAPP™ HINGED KNEE ORTHOSIS - Anterior Closure

PDAC Approved L1832 (ROM Hinged Braces) 633.79/662.17H/530.77UW PDAC Approved L1820 (Standard and Heavy Hinged Braces)

- Constructed of Kuhl[™] perforated neoprene for cool, comfortable compression
- 12" & 16" (30.5cm & 40.6cm) length available in Kuhl neoprene and NoSwet® non-neoprene construction
- 3 hinge styles
- Easily accessible hinge allows interchangeability of hinge styles
- Hinges are fully covered for sports activity
- Elasticized band in thigh strap ensures comfort
- Ideal for pre/post-operative care of MCL/LCL and meniscal repairs
- Longer length can be used for mild to moderate ACL/PCL sprains & strains and as an adjunct to patella fracture

12" Length Knapp **Hinged Knee Orthosis**

Kuhl™ Neoprene

Order No. Description / Options Standard Hinges 5656 5656HH Heavy Hinges 5656ROM Range of Motion

NoSwet®

7656 7656HH 9910.7396.4248 **7656ROM**

Description / Options Standard Hinges Heavy Hinges Range of Motion

16" Length Knapp **Hinged Knee Orthosis**

Kuhl[™] Neoprene

Order No. Description / Options 5658 Standard Hinges 5658HH Heavy Hinges 5658ROM Range of Motion

NoSwet®

Order No. Description / Options Standard Hinges 7658 **7658HH** Heavy Hinges **7658ROM** Range of Motion



Thomas P. Knapp, M.D.

8665.6606.5637

Developed in collaboration with Dr. Knapp, Santa Monica Orthopaedic and Sports Medicine Group, Santa Monica, CA



KNAPP™ HINGED KNEE SIZING:

Measure at mid-patella with leg fully extended

XS: 10" - 12" (25.4cm - 30.5cm) **Med**: 14" - 16" (35.6cm - 40.6cm) **Sm**: 12" - 14" (30.5cm - 35.6cm) 16" - 18" (40.6cm - 45.7cm) 18" - 20" (45.7cm - 50.8cm) **XXL**: 20" - 22" (50.8cm - 55.9cm)

Larger and custom sizes for 5656/5658 models available

upon request



Hinges are removable for easy cleaning



hinge cover allows for easy access



New Range of Motion



Hinge Options:

- 1. Flip Top Pin Stop Range Of Motion Hinge
- 2. Standard Polycentric Hinge
- 3. Heavy Polycentric Hinge







VELOCITY HINGED KNEE

PDAC Approved L1832 633.79/662.17H/530.77UW

- Constructed of Kuhl[™] perforated neoprene to maximize patient comfort and compliance
- Range-of-Motion hinges with stops included; brace can function as a knee immobilizer, a full range of motion brace, or anywhere in between
- Indicated and ideally suited for treatment of ACL, PCL, MCL and LCL sprains, as well as
 post-operative and non-operative meniscus injuries
- Available in anterior closure and pull-up designs
- Hinges are removable for easy brace cleaning
- Can be locked at full extension

Order No. Description / Options

5645 12" (30.5cm) Velocity Hinged Knee - Pull-Up 8555.9895.2727

5646 12" (30.5cm) Velocity Hinged Knee - Anterior Closure

5647 16" (40.6cm) Velocity Hinged Knee - Pull-Up 10695.7847.1909

5648 16" (40.6cm) Velocity Hinged Knee - Anterior Closure





Range-of-Motion Options:

Extension Stops allow 10°, 20°, 30°, 40° and 50° settings, (preset @ 10° upon delivery)

Flexion Stops at 50°, 60°, 70°, 80°, 90° and 100° increments

Brace can be locked at full extension

VELOCITY HINGED KNEE PULL-UP SIZING:


Measure circumference 6" above midpatella with leg fully extended.

XS: 13" - 16" (33.0cm - 40.6cm)
Sm: 15" - 18" (38.1cm - 45.7cm)
Med: 16" - 19" (40.6cm - 48.3cm)
Lg: 18" - 20" (45.7cm - 50.8cm)
XL: 19" - 23" (48.3cm - 58.4cm)
XXL: 21" - 24" (53.3cm - 61.0cm)
XXXL: 22" - 25" (55.9cm - 63.5cm)

VELOCITY HINGED KNEE ANTERIOR CLOSURE SIZING:

Measure circumference 6" above midpatella with leg fully extended.

XS: 12" - 16" (30.5cm - 40.6cm)
Sm: 14" - 18" (35.6cm - 45.7cm)
Med: 16" - 20" (40.6cm - 50.8cm)
Lg: 17" - 22" (43.2cm - 55.9cm)
XL: 19" - 25" (48.3cm - 63.5cm)
XXL: 23" - 28" (58.4cm - 711cm)
XXXL: 25" - 29" (63.5cm - 73.7cm)





5648

5647





Inferior tubular buttress provides patella control

HINGED PATELLA STABILIZER - INFERIOR U TUBULAR BUTTRESS

PDAC Approved L1810 113.43/126.68H/88.14UW

- 12" (30.5cm) Nylon covered neoprene offers comfort, heat and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Neoprene hinge covers allow full sports participation

Order No. Description / Options 5195.7363.6144

3641H Hinged Stabilizer - Inferior U Tubular Buttress, Open Popliteal
Hinged Stabilizer - Inferior U Tubular Buttress, Closed Popliteal

KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 12" - 13" (30.5cm - 33.0cm)
 Sm:
 13" - 14" (33.0cm - 35.6cm)

 Med:
 14" - 15" (35.6cm - 38.1cm)
 Lg:
 15" - 17" (38.1cm - 43.1cm)

 XL:
 16" - 18" (40.6cm - 45.7cm)
 XXL:
 18"-20" (45.7cm - 50.8cm)

Larger and custom sizes are available upon request.



HINGED PATELLA STABILIZER - Universal Tubular Buttress

PDAC Approved L1810 113.43/126.68H/88.14UW

- 12" (30.5cm) Nylon covered neoprene base offers comfort, heat and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Universal tubular buttress is movable/removable for individualized patient care
- Pull up loops aid in application

Order No.

3630H

3730H

Description / Options 4995.7943.6424

Hinged Stabilizer - Universal Tubular Buttress, Open Popliteal

Hinged Stabilizer - Universal Tubular Buttress, Closed Popliteal

SLIP-ON KNEE SIZING: See sizing chart above

- 3630H base with an easy to apply Posterior Closure design
- Constructed of Kuhl[™] perforated neoprene for cool, comfortable compression

Order No. Description / Options 5195.7363.6144

3630HPC Hinged Posterior Closure Stabilizer - Universal Tubular Buttress

POSTERIOR CLOSURE KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 10"-13" (25.4cm - 33.0cm)
 Sm:
 12"-15" (30.5cm - 38.1cm)

 Med:
 14"-17" (35.6cm - 43.1cm)
 Lg:
 16"-19" (40.6cm - 48.3cm)

 XL:
 18"-22" (45.7cm - 55.9cm)
 XXL:
 20"-24" (50.8cm - 61.0cm)

Larger and custom sizes available upon request.



BASIC KNEE SLEEVE

- 1/8" (3mm) Neoprene with 2-sided nylon outer covering
- Provides compression and warmth

Order No. Description / Options

24100 Sleeve - Closed Patella, Closed Popliteal 1285.099.2901

3605 Sleeve - Open Patella, Open Popliteal24105 Sleeve - Open Patella, Closed Popliteal

KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 12"-13" (30.5cm - 33.0cm)
 Sm:
 13" -14" (33.0cm - 35.6cm)

 Med:
 14"-15" (35.6cm - 38.1cm)
 Lg:
 15"-17" (38.1cm - 43.1cm)

 XL:
 16"-18" (40.6cm - 45.7cm)
 XXL:
 18"-20" (45.7cm - 50.8cm)

Larger and custom sizes are available upon request.





KNAPP™ SLEEVE - ANTERIOR CLOSURE

- 1/8" (3mm) Perforated neoprene construction for cool comfortable compression
- Easy application, even over clothes
- Allows extra adjustment for sizing and use with dressings
- Does not require shoe removal to apply

Order No. Description / Options

5607 KUHL[™] Knapp[™] Sleeve 1895.7231.1161

KNEE SIZING:

Measure at mid-patella with leg fully extended.

 XS:
 12"-13" (30.5cm - 33.0cm)
 Sm:
 13" - 14" (33.0cm - 35.6cm)

 Med:
 14"-15" (35.6cm - 38.1cm)
 Lg:
 15" - 17" (38.1cm - 43.1cm)

 XL:
 16"-18" (40.6cm - 45.7cm)
 XXL:
 18" - 20" (45.7cm - 50.8cm)

Larger and custom braces available upon request.



1 205 3 206 4 207

2595.7181.6022 2695.7881.1922 2795.7591.6732 2990.3902.2452

KNEE IMMOBILIZER - *Elastic Straps*PDAC Approved L1830 101.83/169.15H/80.65



- Made from trimable 3/8" (9.5mm) pile foam laminate
- Reinforced, moldable posterior stays provide additional rigidity
- Medial and lateral movable stay pods offer extra support
- Contact closure application straps
- Elastic adjustment straps provide added compression
- Pile exterior is hook & loop sensitive for easy adjustment
- 3/8" (9.5mm) pile covered foam combined with alternating elastic straps make the brace easy to apply and adjust

 Order No.
 Description / Options
 Order No.
 Description / Options

 201 NA
 12" (30.5cm) Knee Immobilizer
 205
 20" (50.8cm) Knee Immobilizer

 203
 16" (40.6cm) Knee Immobilizer
 206
 22" (55.9cm) Knee Immobilizer

 204
 18" (45.7cm) Knee Immobilizer
 207
 24" (61.0cm) Knee Immobilizer

KNEE IMMOBILIZER SIZING:

Fits mid-patella circumference 12"-23" (30.5cm - 58.4cm)



2395.7761.6302 2595.7181.6022 2695.7881.1922 2895.7202.1642 3095.7612.1362 3395.7732.6882

THREE PANEL KNEE IMMOBILIZER PDAC Approved L1830 101.83/169.15H/80.65



- Provides rigid support for the knee
- 3/8" (9.5mm) pile foam laminate with black flannel interior ensures patient comfort
- Two detachable panels for circumferential adjustments for custom fitting
- Posterior and medial/lateral stavs are malleable
- Optional patella strap provides extra knee control
- Available in 6 lengths

Order No.Description / OptionsOrder No.Description / Options21112" (30.5cm) Knee Immobilizer21520" (50.8cm) Knee Immobilizer21316" (40.6cm) Knee Immobilizer21622" (55.9cm) Knee Immobilizer21418" (45.7cm) Knee Immobilizer21724" (61.0cm) Knee Immobilizer

KNEE IMMOBILIZER SIZING:

Fits mid-patella circumference 12" - 24" (30.5cm - 61.0cm)



WEBLY® ZAP™ ANKLE ORTHOSIS PDAC Approved L1902 79.86/88.25H/69.69UW

- Limits ankle inversion/eversion while not restricting plantar and dorsiflexion
- Figure-eight stabilizing straps replicate ankle taping
- Low profile brace fits in shoes and work boots
- Useful prophylactically and for ankle sprains
- Seamless design maximizes comfort
- Elastic strap closure improves compression and gives a clean appearance
- Fits right or left

Order No. Description / Options

314 Webly® ZAP™ Ankle Orthosis 1995.7931.6961

Available in black and white

WEBLY ZAP SIZING:

Men	6-7	7.5-9	9.5-11	11.5-13	13.5-15	15+
Narrow	XXS	XS	S	M	L	XL
Regular	XS	S	M	L	XL	XXL
Wide	S	M	L	XL	XXL	XXL
Women	5-6	6.5-8	8.5-10	10.5-12	12.5-14	
Narrow	XXS	XS	S	M	L	
Regular	XS	S	M	M	L	
Wide	S	M	L	L	XL	



For small ankles and/or heels move down one size. For large ankles and/or heels move up one size. 308 stays may be used for additional support. Stays may be stacked upon one another for additional support.



WEBLY® ANKLE ORTHOSIS

PDAC Approved L1902 79.86/88.25H/69.69UW

- Constructed of heavy gauge ballistic nylon, with non-stretch, non-fraying nylon straps
- Limits inversion/eversion while not restricting plantar and dorsiflexion
- Neoprene tongue and Achilles relief padding aid in player/patient comfort
- Straps sewn to eliminate Achilles irritation
- Useful prophylactically and for chronic ankle sprains
- Low profile brace fits in shoes and work boots

Order No. Description / Options

304 Webly® Ankle Orthosis 1995.7931.6961
 308A Optional Stays (pair) 295.702.152

WEBLY ANKLE SIZING:

	Women	's Shoe Size - US (Europe)	Men's Sho	e Size - US (Europe)
XXS:	5 - 6	(35 - 36)		
XS:	6 - 8	(36 - 38.5)	6 - 7	(38.5 - 40)
Sm:	8 - 10	(38.5 - 41)	7 1/2 - 9	(40.5 - 42.5)
Med:	10 - 12	(41 - 43)	9 1/2 - 11	(43 - 45)
Lg:	12 - 14	(43 - 46.5)	11 1/2 - 13	(45.5 - 47)
XĽ:			13 1/2 - 15	(47.5 - 49)
XXL:			14 - 16	(48 - 50)





PDAC Approved L4350 104.06/157.56H/78.04UW



 Promotes edema reduction and supportive protection for acute ankle sprains

Universal shell design - one size comfortably fits left and right

■ Two Air Cell Design aids in active edema reduction

Heel pad widths can easily be adjusted externally

Welded Straps and Reinforced Rivet ensure strap durability

Thicker heel pad for durability and support

Air cell pressure is adjustable for a better fit

Order No. Description / Options

307 Air Ankle Stirrup 2095.7641.1871

AIR ANKLE STIRRUP SIZING:

Universal

Adjust heel strap for comfort - Fits up to 3" (8cm) width



PSO - PLANTAR STRETCHING ORTHOSIS

PDAC Approved L4396 150.37/155.55UW

- Indicated and ideally suited for treating Plantar Fasciitis, plantar flexion contractures,
 Achilles Tendonitis, heel spurs and other lower extremity overuse injuries
- Comfortable foam padded foot and calf straps are adjustable with easy close buckles
- Ventilated calf area for patient comfort and greater breathability
- Lightweight design and low profile injection molded shell assist in patient comfort and durability
- Soft, removable and washable foam liner provides a comfortable cushion
- Slip-resistant tread (not recommended to ambulate in this device)
- May be used as a night splint

Order No. Description / Options

333 Plantar Stretching Orthosis 4895.7243.1614

PSO SIZING:

 Women's Shoe Size - US (Europe)
 Men's Shoe Size - US (Europe)

 Sm:
 5 - 7 (35 - 37.5)
 4 - 6 (35.5 - 39)

 Med:
 7 - 10 (37.5 - 41)
 6 - 9 (39 - 42.5)

 Lg:
 10+ (41+)
 9+ (42.5+)



PDAC Approved L3170 46.39/48.04UW



- The heel cups are easy to clean, durable and comfortable to wear
- Heel Spur Pad offers added pressure relief

Order No. Description / Options

327 Silicone Heel Cups 1995.7931.6961



SILICONE HEEL CUP SIZING:

Sm:Fits most children's shoesMed:Fits most women's shoesLg:Fits most men's shoes